

PRODUCT FOCUS:

BED OCCUPANCY SENSOR



TASK Community Care | Product Focus

Bed Occupancy Sensor



Part of the TASK CareLink Falls and Memory Support Packages, for service users who have a substantial or critical need. Please contact us for further details on this or other TASK CareLink Call Packages

For more information on the Bed Occupancy Sensor contact us at:

TASK Community Care
Unit 1, Block 3
CityNorth Business Park
Stamullen
County Meath
IRELAND

info@taskltd.ie
Telephone: +353 (0)1 8435889
www.taskltd.ie

TASK Community Care
The Innovation Centre
Northern Ireland Science Park
Queen's Road
BELFAST
BT3 9DT

info@taskltd.com
Telephone: 0845 304 5535
www.taskltd.com

Overview

The Bed Occupancy Sensor is easy to install as it is simply placed under the mattress of a bed. It can be programmed to send an alert immediately when an individual exits the bed. Alternatively, the alarm can be set to send an alert only when a person does not return to bed within a reasonable time.

This means an onsite carer or staff in our 24hr-monitoring centre can take appropriate action to ensure the individual is safe. It can be used as a single device or installed as part of an integrated telecare system.

How it works

The Bed Occupancy Sensor can operate in a number of ways specific to an individual's unique set of requirements.

1. Firstly, the system can automatically sound an attached buzzer as soon as a person leaves the bed. An onsite carer can then ensure the individual is safe.
2. Similarly, it can be used with a radio trigger. This allows an alert to be sent up to 200meters in-house or to a neighbouring home.
3. Finally, with any of the above options or as an alternative, the system can be used with one of our CareLink two-way speech units. This means staff at our 24/7 CareLink monitoring centre will ensure that help is called for the individual if needed.

An additional option is to automatically switch on a bedside light when getting up during the night and switch back off when an individual returns to bed.

Alternatively, if you have an existing NurseCall system in place we can supply fully compatible replacement Bed Occupancy Sensors.

Benefits

Reduces the possibility of an individual falling and being left lying on the floor indefinitely.

Warranty

Bed Occupancy Sensor Pad – 1 Year Warranty
CareLink two-way speech unit – Lifetime Warranty

© Copyright TASK Community Care 2014. Product catalogue printed 07.2014. Please note: Actual product may differ slightly from image shown.